



## UPS Healthcare

Quality Focused. Patient Driven.

### UPS Laboratory Paks

Leak-resistant paks designed with tough material, a tight seal and an absorbent pad that can hold up to 350ml of liquid.

### UPS Laboratory Boxes

Easy-to-assemble boxes available in three distinct sizes so you can find the right fit to help protect your specimens in transit.

# UPS® Laboratory Packaging

## Solutions for Specimen Logistics



**UPS® Laboratory Paks and UPS® Laboratory Boxes** are part of a vast portfolio of UPS healthcare logistics solutions designed to help your specimen shipments arrive at the right location, at the right time, and in the right condition.

- For overnight shipments, UPS Laboratory Paks and Laboratory Boxes help to provide extra layers of protection.
- The one-time use packages are durable and designed to be highly visible.
- This packaging is specifically designed as overwraps for triple-packed Biological Substance, Category B (UN3373) and Exempt Human and Animal Substances.
- This high-visibility, durable packaging helps protect specimens in transit.

Every specimen is important.

# Solutions for Specimen Logistics

## Exempt Human or Animal Specimen

| PACKAGING          | SIZE   | DIMENSIONS     | ORDER NUMBER |
|--------------------|--------|----------------|--------------|
| UPS Laboratory Pak | Single | 18" x 14"      | 01019521401  |
| UPS Laboratory Box | Small  | 13" x 11" x 2" | 010196120    |
|                    | Medium | 16" x 11" x 3" | 010196140    |
|                    | Large  | 18" x 13" x 3" | 010196150    |

## Biological Substance, Category B

| PACKAGING          | SIZE   | DIMENSIONS     | ORDER NUMBER |
|--------------------|--------|----------------|--------------|
| UPS Laboratory Pak | Single | 18" x 14"      | 01019521501  |
| UPS Laboratory Box | Small  | 13" x 11" x 2" | 010196121    |
|                    | Medium | 16" x 11" x 3" | 010196141    |
|                    | Large  | 18" x 13" x 3" | 010196151    |

## Medical Envelope

| PACKAGING        | ORDER NUMBER |
|------------------|--------------|
| Medical Envelope | 010196122    |



## How to Order

UPS Laboratory packaging can be ordered at no cost on [ups.com](https://www.ups.com) using your shipping account number or by contacting your UPS account representative.

## How to Ship

### DROP-OFF

Shipments of both Exempt Human or Animal Specimens, as well as UN3373 - Category B are both accepted at: The UPS Store®, UPS Customer Counters and UPS Drop Box Locations.

### SCHEDULE A PICKUP

To schedule a pickup, call our U.S. Lab support line at 1-844-4UPS-LAB (1-844-487-7522).

### LEAVE AT YOUR UPS LAB SPECIMEN PICKUP POINT

For select locations, UPS works with central labs to set up UPS Lab Specimen Pickup Points in hospitals, clinics, and medical offices. With a designated Pickup Point, you can get consistent pickups of your specimen shipments at a convenient location in your office.



## To Learn More

If you have any questions about shipping biological substances with UPS, please call the UPS Hazardous Materials Support Center at **1-800-554-9964** or visit [ups.com/healthcare](https://www.ups.com/healthcare).